

SPECIFICATION

Product Name	Vibrating Massage Ball
Product Model	ZK-MSG-10933
Material	Silicone Rubber
Speed	2000 RPM
4 Speed Setting	Low, Medium, High, and Alternate Vibration
Power Source	2000mAh Rechargeable Battery
Output (Required but not included)	5V 2A USB power adapter
Included Components	Vibrating Massage Ball, USB Type-C Cable
Product Weight	1.14 lbs 0.52 kg
Product Diameter	6.6 in (w) x 3.5 in (h) 16.9 cm (w) x 8.9 cm (h)

SAFETY GUIDELINES

- Usage Instructions: Apply gentle pressure for up to 1-2 minutes per area.
- Target Areas: Use on soft tissue at high speeds; avoid the head, bones, and joints.
- Discontinue Use: Stop immediately if you experience pain or discomfort.
- Care & Maintenance: Keep the device dry, avoid submersion in water, and supervise while charging.

PAGE 3

ZeenKind
Healthcare

VIBRATING MASSAGE BALL



www.zeenkind.com
support@zeenkind.com

HOW TO USE



1. Power On: Press and hold the power button for 3 seconds to turn it on. The device will start at Low Vibration.

2. Adjust Vibration Speed: Short-press the power button to switch the speed from Low to Medium, High or Alternate Vibration.

3. Start Massaging: Gently apply the massager to the desired areas of your body for relief.

4. Power Off: Press and hold the power button for 3 seconds to turn it off.

CHARGING GUIDE



To ensure proper charging, charge the Peanut Massager using the included USB Type-C cable with a **5V 2A adapter** (required but not included).

LIGHT INDICATOR

	Vibration	Charging
	Press power button once: 1st light turns on – Low Vibration	1st light blinking: Low battery
	2nd Press: 2nd light is on – Medium Vibration	1st + 2nd lights blinking: 50% battery capacity
	3rd Press: 3rd light turns on – High Vibration	
	4th Press: All lights on – Alternate Vibration (low-medium-high cycle)	All lights on: Fully charged